



The Ultimate Defense Against Deadly Superbugs

Revolutionizing the Final Step of Infection Prevention

Candida auris is a deadly superbug first identified in 2009,¹ and it has rapidly become an emerging global threat to patient health. It is often highly resistant to multiple major antifungal drug classes,² and the reported mortality rate for *C. auris* infections is between 30% and 60%.³ Following several consecutive years of rapid annual increases in U.S. clinical prevalence (2018-2019: 44% increase, 2019-2020: 59% increase, 2020-2021: 95% increase),⁴ *C. auris* is a grave and growing concern for patients, clinicians, and infection preventionists.

PrevaHex^{CHX}® and competitive dressings were tested for their ability to inhibit growth of *C. auris* in an *in vitro* zone of inhibition study. Dressing samples (10mm x 10mm) were placed directly into petri dishes of growth media which had been pre-inoculated with *C. auris* (ATCC# MYA-5001). Samples were incubated for 1 day before zone size was measured for each article. This study was conducted at an independent, contract laboratory under GLP and each dressing type was tested in quadruplicate.*



You've Done Everything Right, Now Protect the Site

New PrevaHex^{CHX}® antimicrobial dressings incorporate groundbreaking advanced technology utilizing pure chlorhexidine, providing rapid, broad-spectrum protection of wound and IV catheter sites.


PrevaHex^{CHX}® Antimicrobial Dressings:

- are the first & only dressings that utilize pure chlorhexidine (CHX) offering antimicrobial protection throughout the transparent areas of the dressing
- provide rapid onset of antimicrobial action against a broad spectrum of microorganisms, achieving a high level of efficacy (>5.0 log₁₀ microbial reduction) within 24 hours, and maintaining that high level of efficacy for up to 7 days*
- are effective against gram-positive bacteria, gram-negative bacteria, and yeasts commonly found at wound and IV catheter sites
- demonstrate rapid bacterial elimination of multiple drug-resistant bacterial species both underneath and several millimeters beyond the perimeter of the dressing in *in vitro* laboratory studies*
- incorporate CHX directly into the transparent areas of the dressing which provides the dual benefit of both protecting the site and providing full site visualization
- are offered with integrated securement features to stabilize & secure IV catheters
- are available with a 3-year shelf-life in both single sterile format as well as bulk, non-sterile product format (for use in pre-packed, procedural kits)



a member of the entrotech group of companies

Summary of Results from Comparative Dressing Test Against *Candida auris*:



- PrevaHex^{CHX}® demonstrated a zone of inhibition 50% wider than 3M® Tegaderm® Antimicrobial I.V. Advanced CHG
- PrevaHex^{CHX}® provided an area of protection 90% larger than 3M® Tegaderm® Antimicrobial I.V. Advanced CHG
- 3M® Tegaderm® and Covalon® SurgiClear® showed no efficacy

DRESSING TYPE	ANTIMICROBIAL TECHNOLOGY	ZONE OF INHIBITION BEYOND SAMPLE (mm)*	SURFACE AREA OF PROTECTION BEYOND SAMPLE (mm)*
PrevaHex ^{CHX} ® (1485CH)	CHX: 10% wt./wt.	18mm	528 mm ²
3M® Tegaderm® Antimicrobial I.V. Advanced CHG (9132)	CHG: 2% wt./wt.	12mm	278 mm ²
3M® Tegaderm® (1626W)	N/A	None	None
Covalon® SurgiClear® (TWBD1019)	CHA: 3% wt./wt. Silver Salts: 0.5% wt./wt.	None	None

*Note: *In vitro* effectiveness does not predict clinical performance. (Data on File: entrotech life sciences, llc, Columbus, OH)

References:

1. Satoh K, Makimura K, Hasumi Y, Nishiyama Y, Uchida K, Yamaguchi H. *Candida auris* sp. nov., a novel ascomycetous yeast isolated from the external ear canal of an inpatient in a Japanese hospital. *Microbiol Immunol*. 2009 Jan;53(1):41-4. 2. Ostrowsky B, Greenko J, Adams E, Quinn M, O'Brien B, Chaturvedi V, Berkow E, Vallabhaneni S, Forsberg K, Chaturvedi S, Lutterloh E, Blog D; C. auris Investigation Work Group. *Candida auris* Isolates Resistant to Three Classes of Antifungal Medications - New York, 2019. *MMWR Morb Mortal Wkly Rep*. 2020 Jan 10;69(1):6-9. 3. Chowdhary A, Jain K, Chauhan N. *Candida auris* Genetics and Emergence. *Annu Rev Microbiol*. 2023 Sep 15;77:583-602. 4. Lyman M, Forsberg K, Sexton DJ, Chow NA, Lockhart SR, Jackson BR, Chiller T. Worsening Spread of *Candida auris* in the United States, 2019 to 2021. *Ann Intern Med*. 2023 Apr;176(4):489-495.



To learn more about PrevaHex^{CHX}® antimicrobial dressings, please contact us at info@prevahexchx.com or visit prevahexchx.com

entrotech life sciences, llc
1245 Kinnear Rd.
Columbus, OH 43212
1 (415) 513-4494
www.entrotechlifesciences.com

PrevaHex^{CHX}® is a registered trademark of entrotech, inc.; SurgiClear® is a registered trademark of Covalon Technologies Ltd.; Tegaderm® is a registered trademark of Solventum Intellectual Properties Company.

© 2025 entrotech life sciences, llc. All rights reserved.

MK-20272/Rev. A/022525